

BACKGROUND

Tumor Necrosis Factor alpha (TNF- α) is an inflammatory cytokine secreted by macrophages, monocytes, neutrophils, T-cells, NK-cells following their stimulation by bacterial LPS. TNF- α activates signals through two receptors, TNF-R1, which is expressed on most cell types, and TNF-R2, which is expressed mainly on immune cells. TNF- α can have many functions including, to stimulate of phagocytosis in macrophages, to chemoattract neutrophils, to increase insulin resistance and to induce fever.

Recombinant mouse TNF- α is a non-glycosylated protein, containing 156 amino acids and having a molecular mass of 17.3 kDa.

Alternative Names:

TNFSF2, Cachectin, DIF, Necrosin, Cytotoxin, Cachexin, TNF

Amino Acid Sequence:

MLRSSQNSS DKPVAHVAN HQVEEQLEWL SQRANALLAN
GMDLKDNLV VPADGLYLVY SQVLFKQGC PDYVLLTHTV
SRFAISYQEK VNLLSAVKSP CPKDTPEGAE LKPWYEPYIL
GGVFQLEKGD QLSAEVNLPK YLDFAESGQV YFGVIAL

TECHNICAL INFORMATION

Source: *Streptomyces lividans*

Physical Appearance:

Sterile Filtered white lyophilized (freeze-dried) powder.

Formulation:

Each mg of recombinant mouse TNF- α contains 50mM Tris-HCl pH 8 and 20mM NaCl.

Stability:

Lyophilized product is very stable at -20°C. Reconstituted material should be aliquoted and frozen at -20°C. It is recommended that a carrier protein (0.1% HSA or BSA) is added for long term storage.

Reconstitution:

Centrifuge vial before opening. When reconstituting the product, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution. It is recommended to reconstitute the lyophilized product with sterile water at a concentration of 0.1 mg/ml, which can be further diluted into other aqueous solutions.

Protein Content and Purity determined by:

- UV spectroscopy at 280 nm
- RP-HPLC calibrated against a known standard
- Quantitation against a known standard via reducing and non-reducing SDS-PAGE gels.

Endotoxin Level:

Endotoxin level, as measured by LAL analysis, is <0.01ng/ug or <0.1EU/ug.

Biological Activity:

The activity is determined by the cytolysis of mouse L929 cells in the presence of Actinomycin D and is typically less than 0.1 ng/ml.

Products are for research use only. They are not intended for human, animal, or diagnostic applications.

